TEST HELIX
HELIX TEST

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- 1 - Open the device by unscrewing the capsule containing the indicator strips
- 2 - Fold the strip on itself.
- 3 - Insert the strip into the sterilization chamber. The folded side must be turned towards the outside of the sterilization chamber (steam inlet).
- 4 - Screw on the tip, checking that the tip is correctly positioned on the thread.
- 5 - Start a sterilization cycle. The parameters required to obtain the color change are :
 - 134°C saturated steam for 3.5 minutes
 - 121°C saturated steam for 15 minutes
- 6 - When the cycle is finished, open the device chamber and remove the indicator.
- 7 - When the air inside the device has been removed, the saturated steam enters the tube and reaches the chamber. The temperature of 134°C must be maintained for 3.5 minutes (or 121°C for 15 minutes) so that the indicator turns from yellow to black indicating that all parameters have been met.

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CONDITIONS DE STOCKAGE / STORAGE CONDITIONS

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
Shelf life of the test: 5 years from the date of manufacture under proper storage conditions (10 to 30°C and humidity between 30 and 70%)

CONDITIONS DE STERILISATION / STERILIZATION CONDITIONS

Produit non stérile / *Non sterile product.*

VISUEL PRODUIT / VISUAL PRODUCT

DISTRIBUTEUR / DISTRIBUTOR
MEDISTOCK